



More Info

Data sheet



Protective screen front thermop. canteen

The practical hinged and lockable germ guard screen for the thermoport@canteen.

Easy-to-install hinged and lockable germ guard for the thermoport canteen® for the safe serving of food using an accessory set. The accessory set consists of a transparent acrylic glass panel with CNS glass mounts and matching hinges screwed together at the sides. The accessory set is simply screwed to the existing fold-out shelf mounting points with the two mounting plates supplied. By slightly lifting the panel and rotating it 180°, the panel is hinged upwards and locked in place for serving food. The same applies for the reverse process: Lift the panel slightly, rotate it 180° downwards, lock the panel in place. The thermoport® canteen is ready to move again.

For mounting on the centre of the thermoport® canteen on the guest side.

Retrofittable with a CHECK sensor for seamless and automatic temperature documentation during transport. Clear identification and connection to the CHECK CLOUD are possible with the QR code.

TECHNICAL SPECIFICATIONS

Protective screen front thermop. canteen



TECHNICAL SPECIFICATIONS

Dimensions 397 x 722 x 63 mm

Weight 9,2 kg

Order number 85 01 17 18

BENEFITS

Protective screen made of acrylic glass material PMMA, 5 mm, transparent, screwed together with CNS glass mounts.

Hinges and glass mounts at the sides made of chrome nickel steel, material no.: 1.4301 (CNS).

Hinged and lockable germ guard screen.

Easy to install on the existing fold-out shelf mountings for the safe serving of food.

Easy cleaning with detergent diluted in hot water or with hot water - using a high-pressure cleaner.